



SIDE	NAV POST		TRYKK	SKYV	VARME	RISIKO	FOREB. EFFEKT	SÅRGRAD	BEHANDL. EFFEKT	SITTE STABILITET	POSISJONERING
10	4	ROHO® QS HIGH PROFILE®	1 ●●●●●	●●●●●	●●●●●	høy	●●●●●	høy	●●●●●	●●●●●	●●●●●
10	4	ROHO® QS MID PROFILE®	●●●●●	●●●●●	●●●●●	høy	●●●●●	høy	●●●●●	●●●●●	●●●●●
8	3	ROHO® HIGH PROFILE® Single Valve	1 ●●●●●	●●●●●	●●●●●	høy	●●●●●	høy	●●●●●	●●●	●●
8	3	ROHO® MID PROFILE® Single Valve	●●●●●	●●●●●	●●●●●	høy	●●●●●	høy	●●●●●	●●●	●●
16	4 + 5	ROHO® Hybrid Elite® Dual Valve	●●●●●	●●●●●	●●●●●	høy	●●●●●	høy	●●●●●	●●●●●	●●●●●
16	5	ROHO® Hybrid Elite® Single Valve	●●●●●	●●●●●	●●●●●	høy	●●●●●	høy	●●●●●	●●●●●	●●●●●
12	4	ROHO® CONTOUR SELECT®	●●●●●	●●●●●	●●●●●	høy	●●●●●	høy	●●●●●	●●●●●	●●●●●
10	4	ROHO® QS LOW PROFILE®	●●●●●	●●●●●	●●●●●	høy	●●●●●	lav	●●●●●	●●●●●	●●●●●
8	3	ROHO® LOW PROFILE® Single Valve	●●●●●	●●●●●	●●●●●	høy	●●●●●	lav	●●●●●	●●●	●●
24		BASE™ Sequential (vekseltrykk)	●●●●●	●●	●●●	middels	●●●	høy	●●	●●●●●	●●
7	3	ROHO® MOSAIC®	●●●	●●●	●●●●●	middels	●●●●●	høy	●●●	●●	●●
21		ROHO® HARMONY®	●●●	●●	●●●	lav	●●●	ingen		●●	●
18	5	ROHO® AirLITE®	●●	●	●●	lav	●●	ingen		●●●●●	●●
6	2	Etac Pembroke	●●	●	●●●	lav	●	ingen		●●●	●
20	6	ROHO® HARMONY® Overlay	●	●	●	lav	●●	ingen		●●	●
19	6	Etac Trim	1 ●	●	●	lav	●	ingen		●●	●
22		ROHO® MINI-Max®	●	●	●	lav	●	ingen		●	●
23		Etac Herten	●	●	●	lav	●	lav		●	●
26		Posture LITE™	●	●	●	lav	●	ingen		●	●